


BestTM nomos[®] Cyber Sonalis[®] Ultrasound System



The **BestTM Cyber Sonalis[®] Ultrasound Imaging System** provides superior visualization of HDR, LDR, RF or Cryosurgical procedures. Our patented SimulViewTM Technology provides simultaneous “live” views of the prostate in both planes, thereby increasing treatment accuracy and precision.

Additional Features

- ✓ Capable of supporting over 20 different probe configurations, this system can be used in almost every discipline
- ✓ 21.5" touch screen technology
- ✓ Superior HD image resolution for improved procedure accuracy, speed and physician confidence level
- ✓ Patented probe design
- ✓ Sagittal array provides for 140 mm length of view encompassing the bladder, prostate and perineal tissue
- ✓ Advanced modular software design provides for future upgrade path
- ✓ Convertible to desktop or portable without losing any functionality
- ✓ Dual battery backup built-in to protect from power loss



Best Nomos • 1 Best Drive • Pittsburgh, PA 15202 USA
tel: 412 312 6700 800 70 NOMOS fax: 412 312 6701 www.nomos.com

*Certain products shown are not available for sale currently.

A Breakthrough in Ultrasound Technology


Advanced keyboard technology, totally sealed and self-healing antibacterial keyboard with SensoFoil™ Technology

21.5" touch screen with multi-point touch capabilities

Multiple options for Probe Holder positions and gel bottle cups

Electric lift column for a wide range of height options

Rotating support column for easy positioning



Easily accessible ECG and computer I/O ports


Large rear handles for easy moving

Center mounted probe for maximum efficiency

Rear cable management hooks

Power cord input port easily accessible

Independent rotating/locking casters



Sonalis® Transducers We offer probes for the most common types of procedures and imaging, all with the enhanced image quality and high resolution that you would expect from the BEST name in ultrasounds!

8L2A Linear Array



12L5A Linear Array



14L3 Linear Array



15LW4 Linear Array



15LA Linear Array



15L4A Linear Array



16L5 Linear Array



8V3 Phased Array



4V2A Phased Array



5C2A Curved Array



9MC3 Curved Array



8EC4A Endocavity



XY-BI-Plane Phased Array



10EC4A Endocavity



10BP4 Bi-Plane



8BP4 Bi-Plane



8TE3 Trans-esophageal



Pedoff



16HL7 High Frequency Linear Array



Best Nomos • 1 Best Drive • Pittsburgh, PA 15202 USA

tel: 412 312 6700 800 70 NOMOS fax: 412 312 6701 www.nomos.com