



High Mobility and Table Top Film Cassette Holders

The High Mobility Film Cassette Holder is available in two models; one with vertical tilt only, and the other with vertical and horizontal tilt capability. The High Mobility Film Cassette Holder is extremely economical and allows fast, steady positioning. The unit rolls on three low-friction swivel casters with toe-operated brakes.

Vertical and horizontal adjustments are quick and precise. Angular, horizontal and vertical adjustments are secured by a single turn of locking knobs. The holder can accommodate film cassette sizes from 8 in x 10 in through 14 in x 17 in.

The High Mobility Film Cassette Holder features strong metal construction with a baked-on enamel finish that makes it very durable.

Model 30-3003

This High Mobility Film Cassette Holder features vertical tilt only.

Model 30-3004

This High Mobility Film Cassette Holder features vertical and horizontal tilt.

Specifications

Base Dimensions: 61 cm (24 in)
Height: 73 cm to 122 cm (29 in to 48 in)
Reach: 41 cm (16 in)
Weight: 24 kg (52 lbs)

Model FH-201

The tabletop Flexi-Holder 201 is a lightweight, portable film cassette holder that allows cassette positioning without sandbags. It accepts virtually any type of cassette and can hold cassettes below the edge of the x-ray or operating room table. It does not obscure the cassette film area.

Repeatability is accomplished by means of the triangulation scale on the back of the holder. The FH-201 tilts 135°. Its 3/16 in aluminum frame is both lightweight and sturdy.

Operation is simple and fast. Insert the cassette, tighten the knobs and align the angulation scale. Then slide the FH-201 under the patient or tabletop.

Specifications

Base Dimensions: 46.4 cm x 62.25 cm (18.25 in x 24.5 in)
Weight: 3.4 kg (7.5 lbs)