



UPGRADE KIT

THERATRON® EQUINOX™ EXTERNAL BEAM THERAPY SYSTEM

Best® Theratronics, the world's leading supplier of Gamma Teletherapy equipment, is pleased to announce the immediate availability of Theratron® Equinox™ External Beam Therapy System UPGRADE KIT!

Resulting from numerous requests from our valued customers, the Equinox™ Upgrade provides for existing Theratron® 780C, 780E, Elite 80, Elite 100, T1000E, and Phoenix™ units to be upgraded to modern, state-of-the-art technology in a timely and cost-attractive way.

Upgrade Features

- Calculated Arc Speed
- Graphical Control System
- Asymmetric Jaws (optional)
- R&V System Ready (optional)
- Service Log Files
- On-Board Verification
- Motorized Wedge (optional)
- Collision Detection (optional)

Upgrade Includes

- Removing all old controls and electronics, and installing a new control system and covers
- Replacing the old collimator system with the new Equinox collimator
- Replacing the old treatment table with the new Avanza™ Table
- Retaining the head rotation capability

© 2012 Best Medical International

From this ...



... to this!



The Upgrade Kit not only allows for older units to become new again but it will also allow for **additional treatment options** with the following future modalities:

- MLC for Conformal RT (replacement of shielding blocks) (only for 100cm SAD units)
- MLC for "Sliding Window" IMRT (only for 100cm SAD units)

Furthermore, its **Avanza™ Patient Positioning Table** is fully integrated with the Theratron® Equinox™, providing accurate and reproducible patient positioning, along with efficient and comfortable patient set-up.

For more information, please contact our Customer Service Team in Canada at **613-591-2100, ext.1080** or **service@theratronics.ca**

BT/MF 8015 TEUK (2) v.04.2012